



Optical Mineralogy

WS 2012/2013

Course structure

- THEORY OF OPTICS – first 6 or 7 weeks of lectures and practicals
- MINERAL CHEMISTRY – lectures January onwards
- IDENTIFICATION OF ROCK FORMING MINERALS – pracs January onwards
- EXAM – Week beginning February 4th


Website....

www.staff.uni-mainz.de/tjohnson/Internal/optics.html

- Downloads (pdfs) for lectures and practicals
- Mineral identification sheets
- Links to other resources
- Lots more....

There is lots of EXCELLENT freely available material (English and German) on the internet....

Why use a polarizing microscope??



- **Mineral identification**
- Rock identification
- Microstructural/textural investigation

....this leads to (for example)....

- Crystallisation sequence (igneous petrology)
- Deformation history (structural geology)
- Infer reactions & $P-T$ (metamorphism)
- Diagenetic processes (sedimentary petrology)
- Alteration processes (e.g., weathering)
- and it is CHEAP

Minerals you will learn – 20(ish)

Garnet

Olivine

Clinopyroxene

Orthopyroxene

Hornblende

Glaucophane

Biotite

Muscovite

Chlorite

Epidote

Spinel

Cordierite

Andalusite

Sillimanite

Kyanite

Staurolite

Calcite

Plagioclase feldspar

Alkali feldspar

Quartz

What is the nature of light?

- **Particles** or quanta – photons (Newton) ?

...or...


- Electromagnetic **Waves** (Huygens) ?

....it's both....

Particle-Wave duality (e.g. Einstein, de Broglie)

For the purpose of mineral optics, light is best explained as a **WAVE**....

Electromagnetic Radiation




Transverse waves that are mutually perpendicular

E = Electrical field


M = Magnetic field

C = Propagation direction



For crystal optics we only need to consider the electrical field


The Electromagnetic Spectrum



The Visible Spectrum

DISPERSION through a prism shows us that white light consists of a mixture of all colours of the visible spectrum....

Visible light: 750 nm (red) \rightarrow 400 nm (violet)




Colour




Absorption and Colour

- Selective absorption of certain wavelengths → Absorption colour
- The absorption colour is complimentary to the absorbed wavelengths!
- An example: a green mineral (e.g. hornblende):
 - ♦ Red/orange and blue/violet wavelengths are absorbed
 - ♦ Transparent for green light




Note: Very rarely, colour effects are from interference and diffraction

Reflection and Refraction



Reflection and Refraction



1) Reflection:

Angle of Incidence = Angle of Reflection

$$(\theta_{i1} = \theta_{ref1}; \theta_{i2} = \theta_{ref2})$$

2) Refraction:

Angle of Incidence \neq Angle of Refraction (light is 'bent')

$$(\theta_{i3} \neq \theta_{r3})$$

Velocity of light and refractive index

The ratio of the **velocity** of light in a vacuum/air (c) to the **velocity** of light through a material (v) is called the **Index of Refraction** or **Refractive Index** (n):

$$n_m = c/v$$

n is inversely proportional to v

$$(n_{air} = c/c = 1)$$

Light is slowed down when it enters a denser material, so $n_m > 1$. This slowing down reduces the wavelength of light (the energy and frequency stay constant):

$$\text{Velocity of light: } v = f\lambda$$

$$\text{Light energy: } E = hf = hv/\lambda$$

h = Planck's constant; f = frequency (remains constant); λ = wavelength;
 $c \approx 300\,000$ km/s; v = velocity of light in a material ($v < c \rightarrow n > 1$)

Refraction – Snell's Law

$$\sin\theta_i = AB/CB \rightarrow CB = AB/\sin\theta_i$$
$$\sin\theta_r = CD/CB \rightarrow CB = CD/\sin\theta_r$$

$$\rightarrow AB/\sin\theta_i = CD/\sin\theta_r$$

$$\dots\text{but... } v_1 = AB \text{ and } v_2 = CD$$

$$\rightarrow v_1/\sin\theta_i = v_2/\sin\theta_r \rightarrow$$


Snell's Law

$$v_1/v_2 = \sin\theta_i/\sin\theta_r = n_2/n_1$$

$$\dots\text{as } n_2 = c/v_2, \text{ and } n_1 = 1 \dots$$

$$n_2 = v_1/v_2 = \sin\theta_i/\sin\theta_r$$

Note: v_2 is difficult to measure but $\sin\theta_1$ and $\sin\theta_2$ are not



Reflectance


Reflectance, R :

$$R = \frac{I_R}{I_0} = \frac{(n-1)^2}{(n+1)^2}$$

I_0 = Intensity of incident (source) light

I_R = Proportion of reflected light

n = Refractive Index



Lustre

Glassy: $n = 1,3-1,9 \rightarrow R = 1,7-10\%$

Adamantine: $n = 1,9-2,6 \rightarrow R = 10-20\%$

Sub-metallic: $n = 2,6-3 \rightarrow R = 20-25\%$

Metallic : $n > 3 \rightarrow R > 25\%$



Total Internal Reflection and the Critical Angle

Snell's Law: $n_2 \cdot \sin\theta_{\text{int}} = n_1 \cdot \sin\theta_{\text{ext}}$

$$n_2 \sin\theta_{\text{crit}} = n_1 \sin(90^\circ)$$

$$\sin(90^\circ) = 1, n_1 = 1$$


$$n_2 \sin\theta_{\text{crit}} = 1$$

$$\sin\theta_{\text{crit}} = 1/n_2$$

If $\theta > \theta_{\text{crit}} = \text{T.I.R.}$

Example: Diamond

$$n = 2,42 \rightarrow \theta_{\text{crit}} = 24^\circ$$



Dispersion





Dispersion

The refractive index of a material depends on the frequency and/or the wavelength of the light. - i.e. red light, with a longer wavelength, is refracted less than blue light:


$$n = c/v \rightarrow n = c/f\lambda$$

For most minerals, dispersion is small and not a problem.

However, with other minerals it can be a problem as it affects the other optical properties (e.g. sphene, zircon, some amphiboles)....




Polarization



Normal light is NOT polarized and vibrates in all possible directions!

Plane polarized light (PPL) → only vibrates in a single plane. Light can be polarized by:

- ♦ Reflection
- ♦ Double refraction
- ♦ Selective absorption
- ♦ Dispersion




Polarization by reflection

- Reflected polarized light is parallel to the surface
- Maximum polarization occurs when the angle between the directions of the reflected and refracted beams is 90° :

BREWSTER'S ANGLE θ_B

$$\tan\theta_B = n$$






Polarization by selective absorption

Dichroism = preferred absorption of particular oscillation directions:



Tourmaline - a naturally occurring mineral

Polaroid - a type of synthetic plastic sheet used to polarize light




Polarization by double refraction

- In most minerals (all except those of the cubic system), non-polarized light is split into 2 polarized rays
- The rays have different $n \rightarrow \Delta n = \text{BIREFRINGENCE}$
- These rays are mutually perpendicular
- Example: calcite rhomb - light is split into an ordinary ray (o-ray) and an extraordinary ray (e-ray)



The Polarizing Microscope




Objective and Ocular lenses

- Magnification = Objective \times Ocular:
 - ♦ e.g. $40 \times 10 = 400\times$


- Numerical aperture, A:

$$A = n \cdot \sin\theta$$



- Maximum resolution, d:



$$d = \frac{\lambda}{2A}$$




Thin Sections


- Glass slide
- Glue (Epoxy resin)
- Thin rock slice (**$30\text{ }\mu\text{m} = 0,03\text{ mm}$**)
- Glue (**$n = 1,54$**)
- Glass cover slip

Cover slip
Rock slice
Glass slide





The Polarizing Microscope




Orthoscopic Microscopy

Orthoscopic (parallel light) observations can be made in:

- PLANE POLARISED LIGHT (PPL) - with the analyser **OUT**
 - crystal shape/habit
 - colour/pleochroism
 - cleavage/fracture
 - relief, Becke test → refractive index estimation
- CROSSED POLARISED LIGHT (XPL) - with the analyser **IN**
 - birefringence
 - extinction angle
 - twinning and zoning

PPL - Grain size


- Using a ruler, measure the field of view for each objective lens
- This can then be used to measure maximum and minimum grain-size and grain-size ranges....




PPL - Crystal habit (shape)

- Thin sections are 2d cuts through 3d crystals
- Habits dependent on crystal system, the angle of cut and how perfectly formed the crystals are:


★ EUHEDRAL
= idiomorphic



★ SUBHEDRAL
= hypidiomorphic




★ ANHEDRAL
= xenomorphic




PPL - Crystal habit (shape)


★ EUHEDRAL
= idiomorphic



★ SUBHEDRAL
= hypidiomorphic



★ ANHEDRAL
= xenomorphic



Crystal habits

Acicular **Needle-like**

Bladed **Blade-like**

Equant **Length & width roughly equal**








Fibrous **Slender prisms**

Poikiloblastic **With many inclusions**

Prismatic **Elongate, prism-like**



Tabular **Tablet-shaped**

....etc., etc....



PPL - Colour & Pleochroism


- Colour is caused by selective absorption of certain wavelengths
- Colour (absorption colour) must **always** be observed using PPL
- Pleochroism = direction controlled absorption
 - different colours depend on crystallographic orientation measured by rotating the microscope stage



- Plagioclase is colourless
- Hornblende is pleochroic: light green to olive green


Absorption and Colour

- Selective absorption of certain wavelengths → Absorption colour
- The absorption colour is complimentary to the absorbed wavelengths!
- An example: a green mineral (e.g. hornblende):
 - ♦ Red/orange and blue/violet wavelengths are absorbed
 - ♦ Transparent for green light




Note: Very rarely, colour effects are from interference and diffraction

Pleochroic scheme: Biotite





Pale brown N-S



Dark brown E-W

PPL - Cleavage

- How many?
 - e.g., 0, 1, 2
- Angular relationship?
 - e.g., 90° , 60°
- How well developed?
 - Weak, moderate, good




- **Beware** - Fractures can be easy to mistake as cleavage!


PPL - Relief

- The amount that a mineral *stands out* in the section
- Can be absent, low, moderate, high or very high
- Relief is a measure of the **relative** refractive index (Δn) between the mineral and the epoxy (glue)
- Relief can provide an estimate **of n**

Garnet:	$n = 1,72-1,89$
Quartz:	$n = 1,54-1,55$
Epoxy:	$n = 1,54$




Quartz: very low relief




Garnet: high relief

Relief

Relief can be **positive** or **negative**. A mineral can have moderate relief but a refractive index lower than the epoxy (e.g. fluorite):




Relief





Why do we see relief?

Minerals with different refractive indices (n), cause diffraction, refraction and reflection of the light at grain boundaries:




Becke Line



- As you **lower the stage** (i.e. increase the distance between the objective and sample), **the Becke line moves into the mineral of higher relief....OR....**
- **HHH** = Beim **H**erablassen des Tisches wandert die **h**elle Linie in das **h**öherbrechende Mineral.

Becke Line



The Polarizing Microscope

